



Anti-Nhedc2 polyclonal antibody (CABT-BL2618)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Conjugate	Unconjugated
Applications	WB, IF
Buffer	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 50% glycerol.
Preservative	0.09% Sodium Azide
Storage	Upon receipt, store frozen at -20°C.

BACKGROUND

Introduction	Electroneutral exchange of protons for Na(+) and Li(+) across the inner mitochondrial membrane. Contributes to the organellar volume homeostasis. Required for osteoclast differentiation and bone resorption activity.
--------------	---

GENE INFORMATION

Entrez Gene ID	97086
----------------	-----------------------

Protein Refseq [NP_849208.4](#)

UniProt ID [Q5BKR2](#)
