



Anti-NIM1 (aa 144-193) polyclonal antibody (CABT-BL2622)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from a region within internal amino acids 144-193 (SKLHLVMEYA GGGELFGKIS TEGKLSEPES KLIFSQIVSA VKHMHENQII) of Human Serine/threonine-protein kinase NIM1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	167359
Protein Refseq	NP_699192

UniProt ID

[Q8LY84](#)
