



Anti-NIPAL2 (N-terminal) polyclonal antibody (CABT-BL2623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen Synthetic peptide corresponding to a region within N terminal amino acids AAVAPAGPGDSASAALDESLNFTYGAPGAGNGSLSGDWY RRNQIHLFGV) of Human NIPAL2. NP_079035

Isotype IgG

Source/Host Rabbit

Species Reactivity Human

Purification Immunogen affinity purified

Conjugate Unconjugated

Applications WB, ELISA

Cellular Localization Membrane; Multi-pass membrane protein

Format Liquid

Buffer 2% Sucrose, PBS

Preservative None

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [79815](#)

Protein Refseq [NP_079035](#)

UniProt ID [Q9H841](#)
