



Anti-NKAIN1 polyclonal antibody (CABT-BL2626)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human NKAIN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	NKAIN1 is a member of a family of mammalian proteins with similarity to Drosophila Nkain and interacts with the beta;subunit of Na,K-ATPase (ATP1B1; MIM 182330) (Gorokhova et al., 2007 (PubMed 17606467)).
--------------	--

GENE INFORMATION

Entrez Gene ID [79570](#)

Protein Refseq [NP_078798](#)

UniProt ID [Q4KMZ8](#)
