



Anti-NKAPL (aa 179-228) polyclonal antibody (CABT-BL2629)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 179-228 (KTSRSRNKKK RKNKSSKRKH RKYSDSDSNS ESDTNSDSDD DKKRVKAKKK) of Human NKAPL.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	222698
Protein Refseq	NP_001007532

UniProt ID

[Q5M9Q1](#)
