



Anti-NKX1-2 (aa 101-150) polyclonal antibody (CABT-BL2630)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 101-150 (AHQPERWQGV HEGSPEARAV AVGTEESGAE GLPASPGSPG SPRPRRRRAE) of Mouse Nkx1-2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [20231](#)

Protein Refseq [NP_033149](#)

UniProt ID [P42580](#)
