



# Anti-NOL12 (aa 1-161) polyclonal antibody (CABT-BL2636)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 161 of human NOL12 (Q9UGY1)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF
<b>Cellular Localization</b>	Nucleus; nucleolus.
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	PBS, 1% BSA, 20% Glycerol pH: 7.00
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Nucleolar protein 12 is a protein that in humans is encoded by the NOL12 gene.
---------------------	--

# GENE INFORMATION

Entrez Gene ID [79159](#)

---

Protein Refseq [NP\\_077289](#)

---

UniProt ID [A0A024R1M6](#)

---