



# Anti-nompc polyclonal antibody (CABT-BL2639)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                                                                                    |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>          | A synthetic peptide as a part of <i>Xenopus laevis</i> nompC conjugated to an immunogenic carrier protein                                                          |
| <b>Isotype</b>            | IgG                                                                                                                                                                |
| <b>Source/Host</b>        | Rabbit                                                                                                                                                             |
| <b>Species Reactivity</b> | <i>Xenopus laevis</i>                                                                                                                                              |
| <b>Purification</b>       | Whole antiserum                                                                                                                                                    |
| <b>Conjugate</b>          | Unconjugated                                                                                                                                                       |
| <b>Applications</b>       | WB, IHC-P, IHC-Fr                                                                                                                                                  |
| <b>Reconstitution</b>     | Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.                                                                              |
| <b>Format</b>             | Lyophilised                                                                                                                                                        |
| <b>Buffer</b>             | Whole serum                                                                                                                                                        |
| <b>Preservative</b>       | None                                                                                                                                                               |
| <b>Storage</b>            | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles. |

## GENE INFORMATION

Entrez Gene ID [734212](#)

Protein Refseq [NP\\_001089174](#)

---

UniProt ID [Q4A520](#)

---