



# Anti-NOTUM (full length) polyclonal antibody (CABT-BL2641)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant full length protein (Mouse)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Conjugate	Unconjugated
Applications	IHC-P, WB
Cellular Localization	Secreted.
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	NOTUM (notum pectinacetyltransferase homolog (Drosophila)) is a protein-coding gene. Diseases associated with NOTUM include scabies, and hepatocellular carcinoma. GO annotations related to this gene include hydrolase activity.
--------------	--

# GENE INFORMATION

Entrez Gene ID [77583](#)

---

Protein Refseq [NP\\_780472](#)

---

UniProt ID [Q8R116](#)

---