



Anti-NOXRED1 (aa 179-228) polyclonal antibody (CABT-BL2645)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 179-228 (KLLLNHTNIL RPQYQYDEDS VSVWGANKGV IAALQDPTIL QATCPYSPAG) of Human C14orf148.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	122945
Protein Refseq	NP_001106946

UniProt ID

[Q6NXP6](#)
