



Anti-NPTX1 (aa 186-235) polyclonal antibody (CABT-BL2651)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 186-235 (EGKGGPKNDT EERAKIESAL TSLHQRISL EKGQKDNRPD DKFQLTFPLR) of Mouse Nptx1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic vesicle; secretory vesicle.
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	NPTX1 is a member of the neuronal pentraxin gene family. Neuronal pentraxin 1 is similar to
---------------------	---

the rat NP1 gene which encodes a binding protein for the snake venom toxin taipoxin. Human NPTX1 mRNA is exclusively localized to the nervous system.

GENE INFORMATION

Entrez Gene ID	18164
----------------	-----------------------

Protein Refseq	NP_032756
----------------	---------------------------

UniProt ID	A2ACL9
------------	------------------------
