



Anti-NPTX1 polyclonal antibody (CABT-BL2652)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Cytoplasmic vesicle; secretory vesicle.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	NPTX1 is a member of the neuronal pentraxin gene family. Neuronal pentraxin 1 is similar to the rat NP1 gene which encodes a binding protein for the snake venom toxin taipoxin. Human NPTX1 mRNA is exclusively localized to the nervous system.
--------------	---

GENE INFORMATION

Entrez Gene ID [4884](#)

Protein Refseq [NP_002513](#)

UniProt ID [Q15818](#)
