



Anti-NPTX1 (C-terminal) polyclonal antibody (CABT-BL2653)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to C-terminus of human NPTX1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	NPTX1 is a member of the neuronal pentraxin gene family. Neuronal pentraxin 1 is similar to the rat NP1 gene which encodes a binding protein for the snake venom toxin taipoxin. Human NPTX1 mRNA is exclusively localized to the nervous system.
--------------	---

GENE INFORMATION

Entrez Gene ID [4884](#)

Protein Refseq [NP_002513](#)

UniProt ID [Q15818](#)
