



Anti-NPTX1 polyclonal antibody (CABT-BL2654)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Buffer	PBS with 0.09% (W/V) sodium azide
Preservative	0.09% Sodium Azide
Storage	Maintain refrigerated at 2-8 °C for up to 6 months. For long term storage store at -20 °C in small aliquots to prevent freeze-thaw cycles.

BACKGROUND

Introduction	NPTX1 is a member of the neuronal pentraxin gene family. Neuronal pentraxin 1 is similar to the rat NP1 gene which encodes a binding protein for the snake venom toxin taipoxin. Human NPTX1 mRNA is exclusively localized to the nervous system.
--------------	---

GENE INFORMATION

Entrez Gene ID	4884
Protein Refseq	NP_002513

UniProt ID

[Q15818](#)
