



# Anti-NR1H5 (aa 38-87) polyclonal antibody (CABT-BL2657)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A region within synthetic peptide: QEPPYCQYST AQFPPALQSP SLQSHFNTHG LDPQYSGGSW CGLDARESGF, corresponding to N terminal amino acids 38-87 of Mouse NR1H5
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Cellular Localization</b>	Nuclear
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [381463](#)

Protein Refseq [NP\\_941060](#)

---

UniProt ID [E9Q5A6](#)

---