



Anti-NRARP (aa 36-85) polyclonal antibody (CABT-BL2659)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region of Human NRARP within internal amino acids 36-85 (NMTNCEFNVNSFGPEGQTALHQSVIDGNLELVKLLVKFG ADIRLANRDGW)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	441478
Protein Refseq	NP_001004354
UniProt ID	Q7Z6K4