



Anti-NUPL2 (aa 94-124) polyclonal antibody (CABT-BL2680)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 94-124 of Human NUPL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus, nuclear pore complex. Nucleus membrane; Peripheral membrane protein; Cytoplasmic side. Note: Excluded from the nucleolus.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Nucleoporin-like 2 is a protein that in humans is encoded by the NUPL2 gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID [11097](#)

Protein Refseq [NP_031368](#)

UniProt ID [O15504](#)
