



Anti-NXNL2 (aa 72-121) polyclonal antibody (CABT-BL2683)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 72 - 121 (ADGSCQEMLD FMRELHGAWL ALPFHDPYRQ RSLALLPRLE CSGVILAHCN) of Human NXNL2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Nucleoredoxin-Like 2 is a protein that is encoded by the NXNL2 gene.

GENE INFORMATION

Entrez Gene ID [158046](#)

Protein Refseq [NP_001155097](#)

UniProt ID [Q5VZ03](#)
