



Anti-NXPH1 (internal region) polyclonal antibody (CABT-BL2685)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal of human NXPH1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	May be signaling molecules that resemble neuropeptides and that act by binding to alpha-neurexins and possibly other receptors.
---------------------	---

GENE INFORMATION

Entrez Gene ID	30010
-----------------------	-----------------------

Protein Refseq [NP_689958](#)

UniProt ID [P58417](#)
