



Anti-OBP1F (C-terminal) polyclonal antibody (CABT-BL2691)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide to the C terminal region of Rat Aphrodisin conjugated to an immunogenic carrier protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Secreted.
Format	Liquid
Size	100 µl
Buffer	Whole serum, 1X PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

selectively binds odorants and may be involved in carrying odorants within the nasal epithelium to putative olfactory receptors.

GENE INFORMATION

Entrez Gene ID	192267
----------------	------------------------

Protein Refseq	NP_620258
----------------	---------------------------

UniProt ID	P08937
------------	------------------------
