



Anti-OC90 (aa 230-260) polyclonal antibody (CABT-BL2692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|--|
| Immunogen | Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 230-260 of Human OC90. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB |
| Cellular Localization | Secreted. |
| Format | Liquid |
| Buffer | PBS |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles. |

GENE INFORMATION

Entrez Gene ID [729330](#)

Protein Refseq [NP_001073868](#)

UniProt ID [Q02509](#)
