



Anti-OCM2 polyclonal antibody (CABT-BL2693)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to residues in Human OCM2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Tissue culture supernatant
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Tris glycine, 0.88% Sodium chloride, 0.05% BSA, 40% Glycerol, Tissue culture supernatant pH: 7.40
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction	This gene is similar to the oncomodulin gene, a high-affinity calcium ion-binding protein that belongs to the superfamily of calmodulin proteins, also known as the EF-hand proteins.
---------------------	---

GENE INFORMATION

Entrez Gene ID [4951](#)

Protein Refseq [NP_006179](#)

UniProt ID [P0CE71](#)
