



Anti-ODZ3 polyclonal antibody (CABT-BL2699)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Immunogen | A synthetic peptide corresponding to human ODZ3. |
| Source/Host | Goat |
| Species Reactivity | Human |
| Purification | Antigen affinity purification |
| Conjugate | Unconjugated |
| Applications | IHC, IF, ELISA |
| Format | Liquid |
| Buffer | In Tris saline, pH8.4(0.5% BSA, 0.02% sodium azide) |
| Preservative | 0.02% Sodium Azide |
| Storage | Store at -20°C. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|---------------------|--|
| Introduction | TENM3 (teneurin transmembrane protein 3) is a protein-coding gene. Diseases associated with TENM3 include coloboma, and microphthalmia. GO annotations related to this gene include protein heterodimerization activity and protein homodimerization activity. An important paralog of this gene is TENM2. |
|---------------------|--|

GENE INFORMATION

| | |
|-----------------------|-----------------------|
| Entrez Gene ID | 55714 |
|-----------------------|-----------------------|

Protein Refseq [NP_001073946](#)

UniProt ID [Q9P273](#)
