



Anti-OLFML1 (aa 36-85) polyclonal antibody (CABT-BL2703)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 36-85 (YQRFRLVLEQG LEKCTQATRA YIQEFQEFSSK NISVMLGRCQ TYTSEYKSAV) of Human OLFML1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [283298](#)

Protein Refseq [NP_940876](#)

UniProt ID [Q5HYE3](#)
