



# Anti-OLFR464 polyclonal antibody (CABT-BL2725)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                                                                                                                                                                   |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>          | A synthetic peptide from mouse olfactory receptor family 31, subfamily 2 (MOR31-2) conjugated to an immunogenic carrier protein was used as the immunogen.                                                                                        |
| <b>Isotype</b>            | IgG                                                                                                                                                                                                                                               |
| <b>Source/Host</b>        | Rabbit                                                                                                                                                                                                                                            |
| <b>Species Reactivity</b> | Mouse                                                                                                                                                                                                                                             |
| <b>Purification</b>       | Whole serum                                                                                                                                                                                                                                       |
| <b>Conjugate</b>          | Unconjugated                                                                                                                                                                                                                                      |
| <b>Applications</b>       | IHC, WB                                                                                                                                                                                                                                           |
| <b>Format</b>             | Lyophilised                                                                                                                                                                                                                                       |
| <b>Preservative</b>       | None                                                                                                                                                                                                                                              |
| <b>Storage</b>            | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |
| <b>Ship</b>               | This item will be shipped to you at ambient temperature in a lyophilised form.                                                                                                                                                                    |

## GENE INFORMATION

|                       |                             |
|-----------------------|-----------------------------|
| <b>Entrez Gene ID</b> | <a href="#">258407</a>      |
| <b>Protein Refseq</b> | <a href="#">NP_666524.2</a> |