



Anti-OLFR521 polyclonal antibody (CABT-BL2733)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from mouse olfactory receptor family 101, subfamily 2 (MOR101-2) conjugated to an immunogenic carrier protein was used as the immunogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole serum
Conjugate	Unconjugated
Applications	IHC, WB
Reconstitution	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Format	Lyophilised
Preservative	None
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Ship	This item will be shipped to you at ambient temperature in a lyophilised form.

GENE INFORMATION

Entrez Gene ID [258353](#)

