



Anti-OLFR620 polyclonal antibody (CABT-BL2752)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from mouse olfactory receptor family 4, subfamily 1 and 2P (MOR4-1 and MOR4-2P) conjugated to an immunogenic carrier protein was used as the antigen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole serum
Conjugate	Unconjugated
Applications	IHC, WB
Reconstitution	Add 50 µl of distilled water for a final concentration of 1 mg/ml.
Format	Lyophilised
Preservative	None
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Ship	This item will be shipped to you at ambient temperature in a lyophilised form.

GENE INFORMATION

Entrez Gene ID [258808](#)

