



Anti-OPN3 polyclonal antibody (CABT-BL2785)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|--|
| Immunogen | A synthetic peptide as a part of mouse Encephalopsin conjugated to immunogenic carrier protein |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Mouse |
| Purification | Whole antiserum |
| Conjugate | Unconjugated |
| Applications | WB, IHC-P, IHC-Fr |
| Cellular Localization | Membrane; Multi-pass membrane protein. |
| Format | Liquid |
| Buffer | PBS, Whole serum |
| Preservative | None |
| Storage | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles. |

BACKGROUND

Introduction May play a role in encephalic photoreception.

GENE INFORMATION

Entrez Gene ID [13603](#)

Protein Refseq [NP_034228](#)

UniProt ID [Q9WUK7](#)
