



Anti-OPN3 (internal region) polyclonal antibody (CABT-BL2786)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthesized peptide derived from an internal region of human Encephalopsin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, ICC/IF
Cellular Localization	Membrane; Multi-pass membrane protein.
Format	Liquid
Size	100 µg
Buffer	50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	May play a role in encephalic photoreception.
--------------	---

GENE INFORMATION

Entrez Gene ID	23596
Protein Refseq	NP_055137
UniProt ID	Q9H1Y3