



# Anti-OPN3 (C-terminal) polyclonal antibody (CABT-BL2787)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	This antibody is specific for the C Terminus Region of the target protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Buffer	20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0
Preservative	None
Storage	Store at -20°C. Avoid freeze-thaw cycles.

## BACKGROUND

Introduction	May play a role in encephalic photoreception.
--------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">23596</a>
----------------	-----------------------

Protein Refseq [NP\\_055137](#)

---

UniProt ID [Q9H1Y3](#)

---