



Anti-OPN3 polyclonal antibody (CABT-BL2788)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from the C-terminal domain of human Enkephalopsin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Preservative	None
Storage	Maintain at -70°C in undiluted aliquots for up to 7 months after date of receipt. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction	May play a role in encephalic photoreception.
--------------	---

GENE INFORMATION

Entrez Gene ID	23596
Protein Refseq	NP_055137
UniProt ID	Q9H1Y3