



Anti-OR10T2 (internal region) polyclonal antibody (CABT-BL2790)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthesized peptide derived from the internal region of human OR10T2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Cell membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction OR10T2 belongs to the G-protein coupled receptor 1 family. It functions as an odorant receptor.

GENE INFORMATION

Entrez Gene ID [128360](#)

Protein Refseq [NP_001004475](#)

UniProt ID [Q8NGX3](#)
