



Anti-OR2M7 (internal region) polyclonal antibody (CABT-BL2808)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Immunogen | Synthetic peptide derived from an internal region of human OR2M7. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB, ELISA, ICC/IF |
| Cellular Localization | Cell membrane; Multi-pass membrane protein. |
| Format | Liquid |
| Buffer | 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4 |
| Preservative | 0.02% Sodium Azide |
| Storage | Store at -20°C. Stable for 12 months at -20°C |

BACKGROUND

| | |
|---------------------|--|
| Introduction | Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone |
|---------------------|--|

receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.

GENE INFORMATION

| | |
|----------------|------------------------|
| Entrez Gene ID | 391196 |
|----------------|------------------------|

| | |
|----------------|------------------------------|
| Protein Refseq | NP_001004691 |
|----------------|------------------------------|

| | |
|------------|------------------------|
| UniProt ID | Q8NG81 |
|------------|------------------------|
