



Anti-OR42A (aa 220-259) polyclonal antibody (CABT-BL2812)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: AHMSIFAERLRRLGTYPYESQEQKYERLVQCIQDHKVILR, corresponding to amino acids 220/259 of Fruit fly (<i>Drosophila melanogaster</i>) Odorant Receptor 42a
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Integral membrane protein
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [35514](#)

Protein Refseq [NP_523622](#)

UniProt ID [Q9V9I2](#)
