



Anti-OR43A (aa 210-229) polyclonal antibody (CABT-BL2813)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: LGQIGGEEQSEEERFQRLAS, corresponding to amino acids 210/229 of <i>Drosophila melanogaster</i> Or43a.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>Drosophila</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Integral membrane protein
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [35644](#)

Protein Refseq [NP_523647](#)

UniProt ID [P81917](#)
