



Anti-OR85D (aa 107-150) polyclonal antibody (CABT-BL2823)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen Vector coding for a partial recombinant fusion protein: NISRQRKRLT QVVSRLLEELH PQGLAQQEPY NIGHHLSGYS RYSK, corresponding to amino acids 107-150 of Fruit fly (*Drosophila melanogaster*) Odorant Receptor 85d

Isotype IgG

Source/Host Mouse

Species Reactivity *Drosophila*

Purification Whole antiserum

Conjugate Unconjugated

Applications WB

Cellular Localization Cell Membrane

Format Liquid

Buffer 50% Glycerol, Whole serum

Preservative None

Storage Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID [41011](#)

Protein Refseq [NP_524281](#)

UniProt ID [Q9VHQ2](#)
