



Anti-OSCP1 (aa 1-50) polyclonal antibody (CABT-BL2835)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50(MSVRTLPLLF LNLGGEMLYI LDQRLRAQNI PGDKARKDEW TEVDRKRVLN) of Human OSCP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Basal cell membrane. Syncytiotrophoblast in placenta.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	127700
Protein Refseq	NP_659484

UniProt ID

[Q8WVF1](#)
