



Anti-OSTN polyclonal antibody (CABT-BL2837)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 100 to the C-terminus of Human Osteocrin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted.
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Appears to modulate osteoblastic differentiation. Could also function as an autocrine and paracrine factor linked to glucose metabolism in skeletal muscle.
---------------------	---

GENE INFORMATION

Entrez Gene ID [344901](#)

Protein Refseq [NP_937827](#)

UniProt ID [P61366](#)
