



Anti-OTOP3 (aa 543-572) polyclonal antibody (CABT-BL2838)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 543-572 of Human OTOP3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [347741](#)

Protein Refseq [NP_839947](#)

UniProt ID [Q7RTS5](#)
