



Anti-OTOR polyclonal antibody (CABT-BL2840)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human OTOR.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene is secreted via the Golgi apparatus and may function in cartilage development and maintenance. A frequent polymorphism in the translation start codon of this gene can abolish translation and may be associated with forms of deafness. This gene is a member of the melanoma-inhibiting activity gene family. In addition, alternate polyA sites exist for this gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID [56914](#)

Protein Refseq [NP_064542](#)

UniProt ID [Q9NRC9](#)
