



# Anti-PAC1 (full length) polyclonal antibody (CABT-BL2846)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Full length protein produced in E.coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

Introduction	Encodes the alpha-3 subunit of 20s proteasome.
--------------	--

## GENE INFORMATION

Entrez Gene ID [821774](#)

---

Protein Refseq [NP\\_188850](#)

---

UniProt ID [O81148](#)

---