



# Anti-PAIN (cytoplasmic domain) polyclonal antibody (CABT-BL2847)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide from 4th cytoplasmic loop of drosophila Pain conjugated to an immunogenic carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Fruit fly
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, IHC-Fr
<b>Reconstitution</b>	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
<b>Cellular Localization</b>	Membrane; Multi-pass membrane protein
<b>Format</b>	Lyophilised
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [37985](#)

---

Protein Refseq [NP\\_611979](#)

---

UniProt ID [Q9W0Y6](#)

---