



Anti-PBA1 (full length) polyclonal antibody (CABT-BL2855)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein produced in E.coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Encodes 20S proteasome subunit PBA1 (PBA1).
--------------	---

GENE INFORMATION

Entrez Gene ID [829257](#)

Protein Refseq [NP_001190879](#)

UniProt ID [Q8LD27](#)
