



# Anti-PBG1 (full length) polyclonal antibody (CABT-BL2856)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Full length protein produced in E.coli.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Cytoplasmic and Nuclear
<b>Format</b>	Liquid
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** 20S proteasome beta subunit PBG1 (PBG1) mRNA, complete cds.

## GENE INFORMATION

Entrez Gene ID [842098](#)

---

Protein Refseq [NP\\_176040](#)

---

UniProt ID [Q7DLR9](#)

---