



Anti-PCF1 (aa 98-197) polyclonal antibody (CABT-BL2862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: LKREKERQQREQEKKLREQEKIAAKKMKELEKLEKERIRL QEQQRRKEERDQKLREKEEAQRLRQEQILNKERQQLKLNN FFTKGVEKRIAPNENFVADK, corresponding to amino acids 98/197 of S. pombe CAF1 large subunit.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Schizosaccharomyces pombe
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [2540554](#)

Protein Refseq [NP_596048](#)

UniProt ID [Q1MTN9](#)
