



Anti-Corn PEP4 polyclonal antibody (CABT-BL2880)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Phospho-enol-pyruvate Carboxylase (Maize leaves)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Corn
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Dot Blot, WB, ELISA, IP
Cellular Localization	Cytoplasmic
Format	Liquid
Size	500 µl
Buffer	0.15M Sodium chloride, 0.02M Potassium phosphate. pH 7.2
Preservative	0.01% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [542393](#)

Protein Refseq [NP_001105438](#)

UniProt ID [P51059](#)
