



Anti-PF08_0023 (aa 111-210) polyclonal antibody (CABT-BL2881)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Vector coding for a partial recombinant fusion protein, corresponding to internal sequence amino acids 111-210 of Plasmodium falciparum Plasmodium falciparum HP14, Accession No PF08_0023.Target sequence used to make the antibody: GKITINIDVV QKLNKIMDKI NIE
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Plasmodium falciparum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [2655473](#)

UniProt ID [Q8IB92](#)