



Anti-PF14_0491 (aa 1-100) polyclonal antibody (CABT-BL2883)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein: MTLEHIRQKK GEHNNEKLTE YRQHRKTVD GRLCAAAFVQ DDQTYTDCTR ATDPNGITGK EWCYVEVQLI GKGNRDWDYC KGVINYDVVR SKARTFFQSK, corresponding to N terminal amino acids 1-100 of Plasmodium falciparum HP10
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Plasmodium falciparum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [812073](#)

UniProt ID [Q8IKW0](#)