



Anti-PF4V1 (aa 40-104) polyclonal antibody (CABT-BL2884)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within amino acids 40-104 (GDLQCLCVKT TSQVRPRHIT SLEVIKAGPH CPTAQLIATL KNGRKICLDL QALLYKKIIK EHLES) of Human PF4V1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted
Format	Liquid
Buffer	10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID [5197](#)

Protein Refseq [NP_002611](#)

UniProt ID [P10720](#)
