



Anti-PfRh4 (aa 746-845) polyclonal antibody (CABT-BL2885)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: TLYQEIKTNLENIEQKIMQNIQQTNNRLKINIKKIFDQIN QKYDDLTKNINQMNDKIGLRQMNRLKGKYEEIKKANLQ DRDIKYIVQNNDANINNNNN, corresponding to amino acids 746/845 of Plasmodium falciparum PFD1150c
Isotype	IgG
Source/Host	Mouse
Species Reactivity	P. falciparum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	812435
UniProt ID	C0H496